

L. Belenguer¹, M. Visquert,¹ M. Juan-Borrás¹, I. Escriche¹, R. Fernández,² J.V. García²,
A. Arnau³, Y. Jiménez³, J.J. Manclús³, A. Montoya³

¹Instituto Universitario de Ingeniería de Alimentos para el Desarrollo (IUIAD),

³Centro de Investigación e Innovación en Bioingeniería (Ci2B),


Universitat Politècnica de València.

²Advanced Wave Sensors, Valencia (Spain)

SUMMARY

- The control of chemical residues is essential in honey's marketing to ensure **consumer safety**.
- Equipments available in the market reach the required limit of detection (LOD) (e.g **Liquid Chromatography -HPLC**) but are not suitable for in-situ implantation in honey packaging industry. The reasons are:
 - **Require highly qualified staff**
 - **High cost**
 - **Time-consuming**
- Other approaches available in the market, although faster and without the need of qualified staff, do not reach the required LODs.
- **Piezoelectric immunosensors** based on High Fundamental Frequency Quartz Crystal Microbalance (HFF-QCM) combined with an antigen-antibody recognition process, are alternative or complementary technologies suitable for chemical residues control in food.
- The analytical performance of AWS A20 and AWS F20 platforms (based on HFF-QCM technology), has been tested for the detection of pesticide residues (DDT) in honey.
- AWS platforms lead to the following **benefits**:
 - **Lower cost**
 - **Less sample/reagent consumption**
 - **Label-free, direct and real time transduction**
- DDT standard curves performed in laboratory buffer provided:
 - **Sensitivity (estimated as the I₅₀ value): 0.62 µg/L**
 - **LOD in the 0.05 - 0.1 µg/L range.**


SENSING INTERFACE AND SIGNAL TRANSDUCTION STRATEGY



EXPERIMENTAL SETUP AND TEST PLATFORMS AWS A20 AND AWS F20

Scheme of AWS A20 and AWS F20 platforms


HFF-QCM sensor chip and flow cell assembly



RESULTS

HFF-QCM Standard Curve for DDT

- Concentration of immobilized **BSA-DDT5 conjugate: 10 µg/ml**
- Concentration of **LIB-DDT5.25 monoclonal antibody: 0.5 µg/ml**
- Inter-assay **regeneration: 0.1 M HCl**



Analytical Parameters

Analytical performance of the HFF-QCM immunosensor for DDT as compared to previously developed **ELISA** and **SPR immunosensor**

Analytical parameters (µg L ⁻¹)	HFF-QCM immunosensor	ELISA ¹	SPR immunosensor ²
I ₅₀	0.62	0.74	1.06
LOD	0.08	0.15	0.03
LOQ	0.16	0.22	0.11
WR	0.16 – 2.37	0.22 – 2.47	0.11 – 10.07

¹Estimate from *Abad et al., 1997. J. Agric. Food Chem. 45, 3694-3702.*

² *Mauriz et al., 2007. Biosens. Bioelectron. 22, 1410-1418.*

REMARKS

- **A HFF-QCM immunosensor for DDT insecticide** has been developed.
- Efficient biosensing interfaces were achieved by covalent functionalization of the sensor gold electrodes with mixed alkane-thiol self-assembled monolayers (mSAMs).
- **High-sensitivity monoclonal antibodies** were used as the bio-recognition element of the target analyte in a competitive immunoassay format.
- The limits of detection (**LOD**) and quantification (**LOQ**) and the working range (**WR**) of the developed HFF-QCM immunosensor for DDT are in **the same order of magnitude** as those previously reported for **ELISA** and **SPR immunosensor**.

ACKNOWLEDGEMENT

This work was supported by Spanish Ministry of Economy and Competitiveness project **AGL2013-48646-R**