

Technology Note

QCMD in Lipid Research

Quartz Crystal Microbalance with Dissipation (QCMD) in Lipid Research

QCMD is a label-free surface-analytical technique based on a quartz resonator excited to oscillate at its resonance frequency on one or more overtones. Resonators can have various coatings: gold (Au), silica (SiO_2), titania (TiO_2), etc. It works in aqueous media or organic solvents and is therefore widely used for studying solid/liquid interfaces.^{1,2} At each overtone, QCMD measures changes in the resonance frequency and energy dissipation due to the processes occurring at the resonator surface (Figure 1). Examples of such processes include formation of a film or changes in the geometrical or physical properties of the film.

The key feature that makes QCMD useful in lipid research is its ability to distinguish between different geometries and topologies of lipidic assemblies at interfaces, for example, homogenous solid-supported bilayers or monolayers vs. adsorbed liposomes (Figure 2),³ or other structures (such as cubosomes)⁴ without relying on fluorescent or deuterated labels but by relying on the combination of the frequency and dissipation.

Lipid-related QCMD work can be grouped into several topics, with a total of more than a thousand publications:

- Studies focusing on the interactions between lipids and surfaces.^{2, 5, 6, 7, 8}
- Studies focusing on the properties of the lipids, such as their phase behavior,^{9,10,11} adsorbed liposome deformation,¹² etc.
- Studies examining interactions between lipids and membrane-binding proteins,^{13, 14, 15} peptides,^{16, 17, 18} or viruses.¹⁹ Particularly interesting is that QCMD offers a way to study clustering of membrane-bound proteins.²⁰
- Studies focusing on the interactions of lipids with polymers²¹ or with nanoparticles.²²

It is worth highlighting the studies where QCMD was combined with other techniques, such as neutron reflectometry (NR).^{3, 23, 24, 25} QCMD and NR provide complementary information but




Figure 1. QCMD.

(a) Mineral quartz. (b) A 14 mm 5 MHz QCMD sensor consisting of a quartz disc with gold electrodes. (c) A side view of the sensor depicting its shear motion. Gold electrodes are shown in yellow, sample film in grey, and liquid above in light blue. Thin waves depict the fundamental, thick ones depict the 3rd overtone. Comparing black and blue lines, one sees how the wavelength increases (frequency decreases) in the (sensor + film) case (blue) as compared to the sensor alone (black). The effect of energy dissipation in the film on the attenuation of the amplitude of the wave is shown in red for the 3rd overtone.




Figure 2. Chemical structure of a lipid (top) with its schematic representation (below), a liposome in solution (middle), and surface-supported lipidic structures with low and high dissipation: lipid bilayers (left) and adsorbed liposomes (right). Note the approximate sizes of the structures.

QCMD is much faster and simpler. It can therefore be used to rapidly screen many conditions while NR is typically performed on select few. Combinations of QCMD with other techniques, such as surface plasmon resonance,²⁶ ellipsometry,¹⁷ atomic force microscopy,^{27, 28} and optical reflectometry,²⁹ have also been reported.





Figure 3: Examples of QCMD applications related to lipid work. Top left: lipid-surface interactions; bottom left: lipid phase behavior; right: interactions with nanoparticles, peptides, proteins, polymers, and complex biological samples such as blood or blood plasma.

References

- ¹ Ilya Reviakine, Diethelm Johannsmann, Ralf P. Richter *Hearing what you can't see and visualizing what you hear: interpreting quartz crystal microbalance data from solvated interfaces*. *Anal. Chem.* 83, 8838 - 8848 (2011).
- ² Technology Note: AWSensors QCMD (www.awsensors.com)
- ³ Craig A. Keller and Bengt Kasemo *Surface specific kinetics of lipid vesicle adsorption measured with a quartz crystal microbalance*. *Biophys. J.* 75, 1397–1402 (1998).
- ⁴ Hsin-Hui Shen, Patrick G. Hartley, Michael James, Andrew Nelson, Hugo Defendi and Keith M. McLean *The interaction of cubosomes with supported phospholipid bilayers using neutron reflectometry and QCM-D*. *Soft Matter*, 2011,7, 8041-8049.
- ⁵ Nicole Hain, Marta Gallego, Ilya Reviakine *Unravelling Supported Lipid Bilayer Formation Kinetics: Osmotic Effects*. *Langmuir* 29, 2282-2288 (2013).
- ⁶ Ilya Reviakine, Fernanda F. Rossetti, Alexander N. Morozov and Marcus Textor *Investigating the Properties of Supported Vesicular Layers on Titanium Dioxide by Quartz Crystal Microbalance with Dissipation (QCM-D) Measurements*. *J. Chem. Phys.* 122, 204711 (2005).
- ⁷ Fernanda F. Rossetti, Marta Bally, Roger Michel, Marcus Textor, Ilya Reviakine *Formation and Patterning of Supported Phospholipid Bilayers on the Surface of a Medically Relevant Material*. *Langmuir* 21, 6443-6450 (2005).
- ⁸ Ralf Richter, Anneke Mukhopadhyay, and Alain Brisson *Pathways of Lipid Vesicle Deposition on Solid Surfaces: A Combined QCM-D and AFM Study*. *Biophys J.* 85(5), 3035–3047 (2003).
- ⁹ Astrid Peschel, Arne Langhoff, Eva Uhl, Aruna Dathathreyan, Susanne Haindl, Diethelm Johannsmann, Ilya Reviakine *Lipid Phase Behavior Studied with a Quartz Crystal Microbalance: A Technique for Biophysical Studies with Applications in Screening*. *J. Chem. Phys.* 145, 204904 (2016). DOI: 10.1063/1.4968215
- ¹⁰ Andreas Wargenau, Nathalie Tufenkji, *Direct Detection of the Gel–Fluid Phase Transition of a Single Supported Phospholipid Bilayer Using Quartz Crystal Microbalance with Dissipation Monitoring*. *Anal. Chem.* 86(16), 8017 (2014).
- ¹¹ P. Losada-Pérez, K. L. Jiménez-Monroy, B. van Grinsven, J. Leys, S. D. Janssens, M. Peeters, C. Glorieux, J. Thoen, K. Haenen, W. De Ceuninck, and P. Wagner *Phase transitions in lipid vesicles detected by a complementary set of methods: heat-transfer*



Technology Note

QCMD in Lipid Research

- measurements, adiabatic scanning calorimetry, and dissipation-mode quartz crystal microbalance.* Phys. Status Solidi A 211(6), 1377 (2014).
- ¹² Ilya Reviakine, Marta Gallego, Diethelm Johannsmann, Edurne Tellechea *Deformation of Liposomes Adsorbing on TiO₂ Studied with Quartz Crystal Microbalance.* J. Chem. Phys. 136, 084702 - 084705 (2012).
- ¹³ Edward S. Parsons, George J. Stanley, Alice L.B. Pyne, Adrian W. Hodel, Adrian P. Nievergelt, Anaïs Menny, Alexander R. Yon, Ashlea Rowley, Ralf P. Richter, Georg E. Fantner, Doryen Bubeck, Bart W. Hoogenboom *Single-molecule kinetics of pore assembly by the membrane attack complex.* Nature Communications 10, 2066, 1 – 10 (2019).
- ¹⁴ Ileana F. Márquez, Pablo Mateos-Gil, Jae Yen Shin, Rosalba Lagos, Octavio Monasterio, Marisela Vélez, *Mutations on FtsZ lateral helix H3 that disrupt cell viability hamper reorganization of polymers on lipid surfaces,* Biochim. Biophys. Acta Biomembranes, 1859(10), 1815-1827 (2017).
- ¹⁵ Nicole Heitzig, Alexander Kühnl, David Grill, Katharina Ludewig, Sebastian Schloer, Hans-Joachim Galla, Thomas Grewal, Volker Gerke, Ursula Rescher, *Cooperative binding promotes demand-driven recruitment of AnxA8 to cholesterol-containing membranes,* Biochim. Biophys. Acta Mol. Cell Biol. Lipids 1863(4), 349-358 (2018).
- ¹⁶ Kathleen F. Wang, Nagarajan, Ramanathan, Terri A. Camesano *Antimicrobial peptide alamethicin insertion into lipid bilayer: a QCM-D exploration.* Colloids Surf B: Biointerfaces, 116, 472-81 (2014).
- ¹⁷ Sara Pandidan, Adam Mechler *Nano-viscosimetry analysis of the membrane disrupting action of the bee venom peptide melittin.* Scientific Reports 9, 10841 (2019).
- ¹⁸ Josefine Eilsø Nielsen, Tania Kjellerup Lind, Abdullah Lone, Yuri Gerelli, Paul Robert Hansen, Håvard Jenssen, Marité Cárdenas, Reidar Lund *A biophysical study of the interactions between the antimicrobial peptide indolicidin and lipid model systems.* BBA Biomembranes 1861(7), 1355-1364 (2019).
- ¹⁹ Parveen N., Rimkute I., Block S., Rydell G.E., Midtvedt D., Larson G., Hytönen V.P., Zhdanov V.P., Lundgren A., Höök F. *Membrane Deformation Induces Clustering of Norovirus Bound to Glycosphingolipids in a Supported Cell-Membrane Mimic.* J. Phys. Chem. Lett. 9(9), 2278-2284 (2018).
- ²⁰ Ixaskun Carton, Alain R. Brisson, Ralf P. Richter *Label-Free Detection of Clustering of Membrane-Bound Proteins.* Anal. Chem. 82(22), 9275-9281 (2010).
- ²¹ Fernanda F. Rossetti, Ilya Reviakine, Gabor Csucs, Fabiano Assi, Janos Vörös, Marcus Textor *Interactions of Poly(L-lysine)-g-poly(ethylene glycol) with Supported Phospholipid Bilayers.* Biophys. J. 87, 1711-1721 (2004).
- ²² Desirè Di Silveio, Marco Maccarini, Roger Parker, Alan Mackie, Giovanna Fragnetod, Francesca Baldelli Bombelli *The effect of the protein corona on the interaction between nanoparticles and lipid bilayers.* J. Coll. Interface Sci., 504, 741-750 (2017).
- ²³ Hanna P. Wacklin *Composition and asymmetry in supported membranes formed by vesicle fusion.* Langmuir 27(12), 7698-7707 (2011).
- ²⁴ Alexandros Koutsioubas, Marie-Sousai Appavou, Didier Lairez *Time-Resolved Neutron Reflectivity during Supported Membrane Formation by Vesicle Fusion.* Langmuir 33(40), 10598-10605, 2017.
- ²⁵ Tania Kjellerup Lind, Hanna Wacklin, Juergen Schiller, Martine Moulin, Michael Haertlein, Thomas Gunther Pomorski, Marite Cardenas, *Formation and Characterization of Supported Lipid Bilayers Composed of Hydrogenated and Deuterated Escherichia coli Lipids.* PLoS One. 10(12), e0144671 (2015).
- ²⁶ Erik Reimhult, Charlotte Larsson, Bengt Kasemo, Fredrik Höök *Simultaneous Surface Plasmon Resonance and Quartz Crystal Microbalance with Dissipation Monitoring Measurements of Biomolecular Adsorption Events Involving Structural Transformations and Variations in Coupled Water.* Anal. Chem. 76(24), 7211-7220, (2004).
- ²⁷ Ralf P. Richter, Alain Brisson *QCM-D on Mica for Parallel QCM-DAFM Studies.* Langmuir 20(11) 4609-4613 (2004).
- ²⁸ Elena Rojas, Marta Gallego, Ilya Reviakine *Effect of sample heterogeneity on the interpretation of QCM(-D) data: impurity effects.* Anal. Chem. 80, 8982-8990 (2008).
- ²⁹ Malin Edvardsson, Sofia Svedhem, Guoliang Wang, Ralf Richter, Michael Rodahl, Bengt Kasemo *QCM-D and Reflectometry Instrument: Applications to Supported Lipid Structures and Their Biomolecular Interactions.* Anal. Chem. 81(1), 349-361 (2009).